

exclusive formula XCELENS © 0,3 % UDOCAINE

The experience and reliability already internationally recognized of research and development swiss team **Xcelens** is proud to introduce an innovative **Filler** that combines an excellent performance and results absolutely "No pain".


Xcelens Extra 3 and Xcelens Extra 4 are formulated with a concentration of 23mg/ml of hyaluronic acid crossed-linked with lowest certified BDDE level on the market and with a percentage of lidocaine at 0.3%.
 Xcelens Extra 3 and Xcelens Extra 4 combine the safety of highly purified

hyaluronic acid and cross-linked with BDDE to the high comfort of lidocaine

contained therein.

• MESH TECHNOLOGY

The most of the fillers present on the market usually get a cross-linking so-called "skein structure", instead the **Xcelens** range is made of highly cross-linked hyaluronic acid with a **cross-linking** so called "**mesh**", **Xcelens Extra 3** and **Xcelens Extra 4** have respectively 3 and 4 hyaluronic acid chains of different lengths in addition to "free" hyaluronic acid present into the innovative **Gel Xcelens**.





3

O O ONTENT
 O O ONCENTRATION HA
 O O ONCENTRATION HA
 DEGREE OF CROSSLINKING
 DEGREE OF CROSSLINKING

TREATMENTS: fine lines, superficial wrinkles, glabellar lines, perioral wrinkles, crow's feet, labial commissure, moderate corrections to the nasolabial folds and lip contour.




4

CONTENT

1ML
23MG/ML
23MG/ML
23MG/ML
23MG/ML
23MG/ML
23MG/ML
23MG/ML
24 24 27G
25 2 4 27G
26 2 6 6 6 LIDOCAINE
23MG/ML

TREATMENTS: Deep wrinkles and lines, jaw contour, chin and cheeks area, nasolabial folds, treatment of depressions, increase of lip volume and contour.





EXTRA 3 is indicated for the treatment of fine lines, e specially for perioral wrinkles, crow's feet, glabellar lines, to penetrate the surface and middle stratum of dermis redefining their Appearance.




EXTRA 4 developed to fill deep and middle wrinkles, chin and cheeks area, jaw contour and small depressed areas. Viscosity, softness, elasticity are the main characteristics to keep a NATURAL appearance after the injection.




EXTRA 3 and EXTRA 4 are often involved in the same areas. EXTRA 3 increases the volume of lips and defines the shape of lip contour instead EXTRA 4 has an volumizing effect for lips acting on the deep line and wrinkles.

immediate incredible RESULTS • DURATION

According to independent research centers in Italy and Switzerland, on a sample of 98 people, aged from 29 months to 72 years, the results obtained on the duration of the treatment allows us to ensure the product's effect has a significantly longer duration than other fillers on the market.




Thanks to the presence of LIDOCAINE 0.3%, the treatment will be safe, efficacious and painless, ensuring maximum comfort for the patient.






XCELENS EXTREME HV® the latest and innovative product with Xcelens exclusive formula to treat the lower third of the face.

XCELENS EXTREME HV® main target is to restore and reshape face contours, including the zygomatic area and more specifically the chin and jawline projection, sculpting and shaping the pogonion furrow (angle G - Sn - Pog).



All comparative pre-clinical tests on face plus clinical studies of **XCELENS EXTREME HV***, highlighted a duration of about 18/24 months; in addition, after 3 months from the implant., 100% of the doctors and 97.7% of the interviewed patients evaluated the face volumetric increase as "perfect" or "much better".



Prospective, randomized and controlled study of 100 patients treated exclusively on the 3rd lower of face from 1 to maximum 3 syringes per treatment (executed in maximum 2 visits - initial and retouch) with an interval of about 90 days. Patients returned for follow-up visits at months 1, 3, 6, 12, 18, 24.



EXTREME MY

XCELENS EXTREME HV° the ultimate available choice to reshape face contours restoring any volumetric deficit, including facial LIPOATROPHY. The exceptional viscoelastic properties of **XCELENS EXTREME HV°** combined with the new ergonomic and glass syringe give to doctors an incredible volumizing product together with an excellent hydrogel extrusion from the syringe.




Patient satisfaction, comfort and durability

•	•	•	• •	CONTENT	1ML
•	•	•	• •	CONCENTRATION HA	25MG/ML
•	•	•	• •	RETICULATING AGENT	BDDE
•	•	•	• •	DEGREE OF CROSSLINKING	*****
				NEEDLES	2 X 27G
•	•		• •	NEEDLES	2 x 2/G

INDICATIONS: maxillary redefinition, zygomatic depressions, chin.

G PRIME [G']: rheology as the absolute judge

Rheology is the science that studies the flow and deformability of matter. G' is a rheological parameter indicating the filler capacity to return to its ordinary shape when exposed to dynamic forces. High value of G'Prime characterizes a hard, robust and compact gel with a better response to compression (greater elasticity). For this reason, a product with a greater G' is most suitable for thicker skin and deeper injections. Also an excellent extrusion of the hydrogel from the syringe.



Values and physical-chemical characteristics, highlighted in the study of main fillers available worldwide, are the result of an independent research performed by Xcelens International.





















follow us

■ Xcelens

© xcelensgenyalofficial

XCELENS INTERNATIONAL Sagi

Via Soldini 22 • Chiasso • CH 6830 • Switzerland ph. +41 786.58.57.10 • info@xcelens.ch • xcelens.ch

