



GENYAL[®]

the ultimate generation of fillers



Xcelens Exclusive Formula

Xcelens International is pleased to present an innovative and exclusive range of hyaluronic acid products, uniquely developed to guarantee medical staff, and therefore patients, total comfort, satisfaction and reliability during filler implantation and bio-stimulation treatments: **GENYAL®**!

The **GENYAL®** range is developed by an expert **Swiss-Italian team** of chemists, technicians, engineers and professionals operating in the field of aesthetic medicine with a view to creating more advanced and technologically avantgarde solutions year after year.

- **Double action biphasic** • **Easy to handle** • **Immediate results** • **Natural effects**
- **Highly purified HA** (Lowest level of endotoxin < 0.05 EU/g) • **Great visco-elasticity balance**
- **Natural reabsorption** • **Long-lasting** (8 to 18 months)

Continual **physico-chemical studies** on cross-linking have resulted in the synthesis of all the experience gained to date, resulting in the creation and development of an **ingenious** product that incorporates all the benefits desired by medical staff and patients, **i.e. safe, natural, long-lasting and comfortable implants**, in a single solution: for this reason Genyal Polyvalent, the polyvalent filler for all types of cosmetic correction, was formulated.

Genyal® Polyvalent:

- treats all types of superficial and medium-depth wrinkles
- remodels the lip contours and increases lip volume.

SPECIAL FEATURES

The latest advances in the field have shown that polyvalent fillers are, currently, the best possible upgrade but Genyal's special features make it a unique product.

- **Biphasic preparation of non-animal origin**
sodium hyaluronate synthesised by bio-fermentation dual-action, high-concentration filler: 23mg/1ml of cross-linked hyaluronic acid with bio-stimulating free hyaluronic acid.
- **Ultra high molecular weight product (3.000.000 Daltons)**
this characteristic, together with cross-linkage, allows us to increase the resistance to degradation by hyaluronidases to a maximum and thus maximise the duration of the product in the tissue.
- **Cross-linking with butanediol diglycidyl ether (BDDE)**
the cross-linking agent universally recognised as the least toxic of all, which ensures that the risk of allergies is minimal.

HYALURONIC DERMAL FILLER
GENYAL®
Polyvalent

the filler for every solution

CONTENT	1ML
CONCENTRATION HA	23MG/ML
RETICULATING AGENT	BDDE
DEGREE OF CROSSLINKING	★★★★
NEEDLES	2 X 27G

TREATMENTS: lip contour, barcode, glabella, facial wrinkles, forehead wrinkles, periorcular wrinkles, naso genian folds, perilabial wrinkles, labial commissures.



Genyal® Volumae is the filler with the best balance between high viscosity, elasticity and long lasting. Combined to total confort, satisfaction and reliability during filler implantation, for an excellent, natural, finished filling effect.

The excellent results achieved with Genyal® Polyvalent led to the creation of a specific product that gives volume to the face while maintaining the characteristics of the Genyal® range: the result is a highly cross-linked gel (also containing 23 mg/ ml of hyaluronic acid) which is injected with a 27G needle to reach the desired level of filling without causing trauma to the tissues.

Genyal® Volumae:

fills lines and treats deep wrinkles in thick skin

restores or creates facial volume

increases the volume of the lips

The duration is exceptional for a reabsorbable product: **it varies between 8 and 18 months following implantation.**

SPECTACULAR EFFECTS AND EXCEPTIONAL DURATION

HYALURONIC DERMAL FILLER

GENYAL®

Volumae

for all types of volume

CONTENT	1ML
CONCENTRATION HA	23MG/ML
RETICULATING AGENT	BDDE
DEGREE OF CROSSLINKING	★★★★★
NEEDLES	2 X 27G

TREATMENTS: cheeks, naso genian folds, labial commissures, increase in lip volume, chin, facial wrinkles, nasal correction.



It is now indisputable that filling wrinkles and increasing the volume in the face is not the only way to obtain total rejuvenation of the face; this is why natural, non-cross-linked hyaluronic acid has an important role in combating and preventing skin aging, allowing skin elasticity and firmness to be regained.

The particular formulation of Genyal® Genyalift, with its extremely high concentration of hyaluronic acid (18 mg/1 ml) together with a high concentration of glycerol (20 mg/ml) and ultra high molecular weight (3 million Daltons: double or even three times that of many rival formulas!), allows the recreation of the optimum physiological environment for fibroblast proliferation and migration and for the neo-synthesis of collagen, elastin and other constituents of the extracellular matrix.

Genyal® Genyalift can be used not only on all areas of the face, but also in other regions such as the neckline and the backs of the hands, with different programmes depending on the degree of skin aging; generally, 2-3 initial treatments are planned at short time intervals (15-30 days) and then other maintenance sessions are continued about 2-3 months apart.

HA-GLYCEROL BIO-RESTRUCTURING GEL

GENYAL®

Genyalift

bio-stimulation at the highest level

CONTENT	1ML
CONCENTRATION HA	18MG/ML
STABILISING AGENT	GLYCEROL 20MG/ML
MOLECULAR WEIGHT	3MLN DALTON
NEEDLES	2 X 30G

TREATMENTS: glabella, forehead wrinkles, periocular wrinkles, cheeks, naso genian folds, perilabial wrinkles, labial commissures, neck, facial wrinkles, hands, stretch marks, décolleté, face.



The new Black

The new **Genyaltix**[®] technology allows us to exploit a new type of cross-linking where the long chains of hyaluronic acid are collected in a helical form in a mesh matrix using a new production process that makes the gel monophasic, therefore totally cross-linked and entirely protected from enzyme attack. Therefore, it does not result in adverse reactions to the hyaluronic acid, such as glaucoma or infections, as long as the doctor respects the sanitary rules and does not inject too close to the surface or accumulate too much product in the same point.

solution for timeless beauty

- **Safety and reliability**

The hydrogel injected into the tissue allows the passage of oxygen by osmosis, thus regularly protecting all the physiological functions of the tissue.

- **High yield and durability**

The **Genyaltix**[®] technology used to crosslink the hyaluronic acid gives the chains contained in the hydrogel strong dynamic resistance: the product injected into the tissues will thus remain aggregated in the tissue and will maintain its shape, minimizing any possibility of migration.

- **Details of monophasic hydrogel**

The unique structure of monophasic hydrogel in the new **Black Label Edition of Genyal**[®] is obtained by the total cross-linking of hyaluronic acid without free parts, creating a product that is very compact and less fluid than the biphasic gel, with a very high degree of cohesiveness that under certain conditions will result in greater resistance to enzymatic action.

Label Edition

by GENYAL[®]




GENYAL[®] 20
HYALURONIC DERMAL FILLER



GENYAL[®] LYPS
HYALURONIC DERMAL FILLER



GENYAL[®] 24
HYALURONIC DERMAL FILLER



GENYAL[®] 26
HYALURONIC DERMAL FILLER



GENYAL[®] 200
HYALURONIC ACID BODY FILLER

Genyal® 20 is the most fluid formulation of the line, a light gel to minimize the pressure on the tissue, with an ideal viscoelasticity for an optimal diffusion and positioning capacity: it reduces the risk of edema, for a natural and homogeneous correction.

Thanks to the innovative Genyal® technology, it rejuvenates, tones, plumps the skin and restores lost volume. Its rheological profile is particularly suitable for the fragile periocular area; the implant is placed to fill the lacrimal sulcus: once filled and increased in volume, the **dark circles** will look fresh, illuminated and rested.

TREATMENTS: fine wrinkles, temple filling, acne scarring, glabellar lines, barcode, corners of the mouth, periorbital area wrinkles, tear through.



Genyal®Lyps fills the glabellar lines, the slightly evident nasolabial furrows, the lines of the cheeks and the superficial labial commissures.

It can also be used for redefining the nose and lips and their profile by restoring volume.

TREATMENTS: glabella wrinkles, peribuccal wrinkles, slightly marked nasolabial furrows, labial commissures, cheek wrinkles, lip contour, rinofilling.



GENYAL® 20

HYALURONIC DERMAL FILLER

CONTENT	1ML
CONCENTRATION HA	20MG/ML
RETICULATING AGENT	BDDE
DEGREE OF CROSSLINKING	◆◆
NEEDLES	2 X 30G

**fine lines
& eyes contour**

GENYAL® LYPS

HYALURONIC DERMAL FILLER

CONTENT	1ML
CONCENTRATION HA	22MG/ML
RETICULATING AGENT	BDDE
DEGREE OF CROSSLINKING	◆◆◆
NEEDLES	2 X 27G

**filling the lips
and the superficial lines**

Genyal® 24 is also effective in plumping out "puppet wrinkles", i.e. pronounced nasolabial furrows, and it can give renewed fullness to the lips, fix cheekbones and to treat all types of wrinkles in general. This very thick product is injected into the medium and deep dermis.

TREATMENTS: deep nasogenic wrinkles, anger lines, in general all types of medium and deep wrinkles.



This particular product exhibits the highest concentration of hyaluronic acid in the line **Genyal** (26 MG / ML) and the finest elastic-plasticity characteristics.

Recommended for increasing volume, redefining the oval and face, facial lipoatrophy, chin and jaw area, with an application indicated therefore in the middle and deep dermis.

TREATMENTS: cheekbones, cheeks and chin, face contour.



GENYAL® 24

HYALURONIC DERMAL FILLER

CONTENT	1ML
CONCENTRATION HA	24MG/ML
RETICULATING AGENT	BDDE
DEGREE OF CROSSLINKING	◆◆◆◆
NEEDLES	2 X 27G

polyhedric and monophasic HA filler

GENYAL® 26

HYALURONIC DERMAL FILLER

CONTENT	1ML
CONCENTRATION HA	26MG/ML
RETICULATING AGENT	BDDE
DEGREE OF CROSSLINKING	◆◆◆◆◆
NEEDLES	2 X 25G

volumizing and face contouring



GENYAL[®] 200
HYALURONIC ACID BODY FILLER


**The filler for a perfect reshaping
of the body**

Genyal® 200 is an innovative, implantable medical device, composed of hyaluronic acid and designed to noticeably re-sculpt the body. With Genyal 200, it is possible to give shape to the buttocks, calves and in general, any area of the body where an increase in adipose tissue is appropriate.

Benefits and improvements: it is now possible to reshape the body in a natural and trustable way with Genyal® 200, is biocompatible and can be reabsorbed. It is long-lasting, thus you can transform the appearance of your body with a process that is fast and temporary, after which you can return to work within a very short time, avoiding the scars and risks of prosthetic implants through surgical operations using that may even use total anaesthesia.

It is recommended to use a 16 or 18G cannula.

TREATMENTS: buttocks, calves, depressions and asymmetries, remodelling body profiles.



GENYAL[®] 200

HYALURONIC ACID BODY FILLER

the body filler

CONTENT	10 ML
CONCENTRATION HA	20 MG/ML
RETICULATING AGENT	BDDE
DEGREE OF CROSSLINKING	● ● ● ●
NEEDLES	-





Exclusive World Distributor
xcelens
international

Xcelens International Sagl
Via Soldini 22, Chiasso • CH-6830
Switzerland • Tel: +41 78 658 57 10
info@xcelens.ch • www.xcelens.ch

follow us:  **Xcelens**  **xcelensgenyalofficial**

